


TCS—200MLA

Electronic Personal Scale With Height Measurement

User Manual



Contents

1、 Main Technical Parameters	2
2、 Assembly Method	2
3、 Method of Application	3
4、 Keys Function Description	4
5、 Select The Operation Of The Unit Of Measurement Method	4
6、 Other Function Operation	5
7、 Announcements	6

TCS- 200MLA series uses -high- precision load cell, the special microcomputer control, high definition large screen LCD display, with functions of accurate weighing, measuring height, use reliable and convenient.

It can be stored for 7 personal's information, included sex, age, body height and BMI, the displayer will show the the health index according to inputting parameters.

The products are widely used in family, medical treatment, health care and other places.

1. Main Technical Performance

1. 1 Specification:


NO.	Model No.	CAPACITY	MIN WEIGHT	GRADUATION	BACKLIGHT
01	TCS- 60MLA	60kg	100g	20g	YES
02	TCS-150MLA	150kg	100g	50g	YES
03	TCS-200MLA	200kg	100g	50g	YES
04	TCS-300MLA	300kg	100g	50g	YES
05	RANG OF MEASUREMENT	MAX:210cm	MIN:87cm	1mm	

1. 2 Two ways of Power supply : AC/DC. AC Input 110-240V, 50HZ and AC output 6V. And with 4V DC rechargeable battery inside of scale.

2. Assembly Method

**2. 1 To stock and ship, this scale is packed by separated into 2 parts. A. scale base
B. displayer and pole (See below drawing)**

- 2.2 Insert the seat of pole and into scale base, insert the pole into the base seat , make the display face inside. Fasten it by screw.
- 2.3 Replace battery: Open the batteries case on bottom of scale, get the battery out, and insert new battery. Cover battery case.



3. STEPS TO USE

3. 1 Put the scale onto a hard, level surface, adjust the footing to make it in level.
3. 2 Shut the switch behind displayer, then turn on scale by short press “ON/TARE+” key, scale self-test. After this, it goes into weigh mode.
3. 3 Stand in the center of the platform, when the number displays, and a circle displays on left top corner , the number is the weight.

3. 4. Pull the high measurement , and make the block touches head, then read the4.

FUCTIONS OF KEYS

4. 1 “ON/TARE+” KEY

A . In off status, Short press this key to turn on scale. B. In weigh mode, short press this key to Zero the weigh to “0”. In setting mode, short press this key to make one increase.

4. 2 “OFF” KEY

Short press this key to turn off scale.

4. 3 “UNIT -” KEY

To set weighing unit, kg or lb. the displayer shows the unit on right side KG and LB. Weighing unit stays in KG when turn on scale.

4. 4 “MEM” KEY

Short press this key to set memory.

4. 5 “SET” KEY

Short press this to insert user information, such as sex, age, high.

4. 6 “RECALL” KEY

Short press this key to display memories.

5. Steps to select weighing unit

5 . 1 Weighing unit can by “selectable” and “locked”. Turn on scale, press “UNIT” key, weighing unit can swift, that is “selectable “, otherwise is “locked”. “Selectable” status can swift weighing unit an any time, “Locked” status remains a weighing unit and can not be changed.

5. 2 Operation to select “selectable” and “locked”: turn on scale, when displays “0”, press and hold “UNIT” key until weighing unit kg or lb flash, the flashing weighing unit is the unit which will be locked. Repeat above steps and swift status between “selectable” and “locked”.
5. 3 Any error operation , please turn off the scale and re-operate above steps.

6. OPERATION TO OTHE FUCTION

6. 1 Memory Function: Turn on scale, weigh and get result, after circle displays on left top corner , short press “MEM” key, displays “SAVE” and memory stocks. This scale can stock 7 users, each users 7 arrange memories. If the user memories more than 7 arrange, the earliest memory will replaced by new memory.
6. 2 Recall Function: Turn on scale, short press “RECALL”, displayer shows “SLUSE”, short press “ON/ZERO+” key can select user, after selected user, short press “RECALL” key, the latest memory will display. Short press “ON/TRADE+” key can display more memories. Short press “RECALL “ key to end recall mode.

Remark: If no any memory was stocked, when recall, it displays “noDAT”

6. 3 Input Information: Turn on scale, short press “SET” key . User number flashes, short press “ON/TARE+” key to increase, “UNIT-“ to reduce, “SET” key to confirm. After this, use same way to set sex, age, high.
6. 4 Clean Up Function: In recall mode, press and hold “ON/TARE+” key until displays ”CLRPX”, this user’s all memories clean up.

REMARK: About BMI

world health organization Standard:

Less than 18.5 thin, 18.5~24.9 normal, More than 25 partial fat,
More than 30 fat, More than 40 Obesity

Asia Standard:

Less than 18.5 thin, 18.5~22.9 normal, More than 23 partial fat,
More than 30 fat, More than 40 Obesity

7. NOTE

7. 1 When works by battery, when displayer shows “LO”, please charge the battery. When charge, please use the wire which can be found in package.
7. 2 Scale should put onto hard, level surface, adjust the footing to make it level to get a precision result .
7. 3 Clean the scale after used, keep it dry.
7. 4 User’s disassembly or un-normal use, (over load, strong impact ect.) leads to load cell broken, is out of guarantee.
7. 5 Any doubt, please contact local retailers or give us phone call.